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- [1] NMIJ-PTB: *Bilateral comparison on nanometrology according to the rules of CCL key comparisons*; APEC ISTWG Project Interlaboratory Comparison on Nanoparticle Size Characterization, 2012
- [2] Projeto NanoValid, *Developing Reference Methods for Nanomaterials*, 2012-2015
- [3] NANoREG, *A common European approach to the regulatory testing on nanomaterials*, 2014-2017. Protocolos disponíveis em: <http://www.nanoreg.eu/>
- [4] C.P.Gouvêia, S.M.Landi, C.E.Galhardo, J.C. Damasceno. Nota técnica Inmetro: *Análise dimensional de nanomateriais utilizando microscopia eletrônica*, em http://www.inmetro.gov.br/metcientifica/pdf/Analise_dimensional_de_nanomateriais_utilizando_microscopia_eletronica_-_Nota_Tecnica.pdf
- [5] ISO 29301:2017. *Microbeam analysis – Analytical transmission electron microscopy – Methods for calibrating image magnification by using reference materials having periodic structure*. 2010, ISO: Geneva, Switzerland.
- [6] J.A.W.M. van der Laak et al.: *Automated magnification calibration in transmission electron microscopy using Fourier analysis of replica images*. Ultramicroscopy 106, pp. 255–260, 2006 [DOI: 10.1016/j.ultramic.2005.08.001]
- [7] NIST: Bonevich, J.E. and W.K. Haller, *Measuring the size of nanoparticles using transmission electron microscopy (TEM)*, in NIST-NCL Joint Assay Protocol, PCC-7 Version 1.1, NIST, Editor. 2010, NIST-NCL: Washington, DC, United States
- [8] NPL: Boyd, R.D., et al., *Good practice guide for the determination of the size distribution of spherical nanoparticle samples*, in Measurement Good Practice Guide No. 119. 2011, National Physical Laboratory: Teddington, Middlesex, UK.
- [9] ISO/TS 80004-2:2015. *Nanotechnologies – Vocabulary – Part 2: Nano-objects*. 2015, ISO: Geneva, Switzerland.
- [10] L 275/38 Official Journal of the European Union 20.10.2011
- [11] Mast, P.-J. De Temmerman. Deliverable D 2.10 *Protocol(s) for size-distribution analysis of primary NM particles in air, powders, and liquids*. CODA-CERVA Veterinary and agrochemical research center. Documento público sob Creative Commons Attribution-NonCommercial-ShareAlike 4.0 International License. 2016.

- [12] P.-J. De Temmerman et al.: *Measurement uncertainties of size, shape, and surface measurements using transmission electron microscopy of near-monodisperse, near-spherical nanoparticles*. J Nanopart Res 16:2177, 2014 [DOI: 10.1007/s11051-013-2177-1]
- [13] H.G. Merkus: *Particle size measurements: fundamentals, practice, quality*. Pijnacker: Springer, 2009. 533 p. ISBN 978-1-4020-9015-8
- [14] ISO 13322-1:2014. *Particle size analysis – Image analysis methods – Part 1: Static image analysis methods*. 2014, ISO: Geneva, Switzerland.

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Elaborado por	Sandra Marcela Landi	Pesquisadora -Tecnologista
Aprovado por	Oleksii Kuznetsov	Chefe da Dimat